



LIST OF MATERIAL			
ITEM	PART NUMBER	QTY.	DESCRIPTION
1	007818	1	HTR TUTCO 1.2KW 120/240V
2	336929	1	MTR MTG LAB LINE .04HP
3	338125	1	HEATER COVER
4	338126	1	LEFT SIDE DUCT
5	338127	1	RIGHT SIDE DUCT
6	336928	1	SLV 2.50 LG ALUM
7	008333	1	MTR .04HP 1575 RPM
8	007281	1	FAN WH 05.75X01.00X0.313 ALZ CCW
9	016643	1	CLIP TINNEMAN C17120-017-4
10	133033	21	SCW SM PHTX ZIN #06X0.50
11	012777	6	WSHR FLT 18-8 #06
12	336959	2	SHIM STRIKE 0.50X0.50 SST
13	008199	2	LATCH #350F2C
14	272516	4	SCW MAC THTX18-8 #08-32X1.00*
15	272579	4	SCW MAC THTX18-8 0.25-20X0.50
16	148048	1	DOOR PULL-AUSTIN #AE1-30357
17	143614	2	SCW TT PHTX ZIN #08-32X0.50
18	272530	4	SCW MAC PHTX ZIN 0.25-20X1.00
19	015229	1	ENCL BLWR UF12A23 STH 230V
20	018967	4	SCW MAC RHSL ZIN #06-32X2.00
21	008484	4	NUT HX ZIN #06-32
22	012814	4	WSHR LK ZIN #06 SHKPRF INT
23	144755	1	T/C JJ38005.5-303.5-25048-006
24			
25			
26			
27			
28			
29			
30	336951	1	2 INCH PORT COVER PLATE
31	336950	1	PORT TUBE 2.0 IN (6.21LG)
32	132623	6	RIV POP SST .12DIAX.032X.187 LG
33	057939	1	STOPPER SLCN #10.5 HI-TEMP*
34	163101	1	T/C K48U-006-00-19-F1B084-6
35	163023	1	POCKET PROBE READOUT
36	272527	2	SCW TT PHTX ZIN #08-32X1.00
37	325729	1	METER DG OMGA HH801B HAND HELD DUAL CHANNEL
38	141894	1	LT BULB 40A15 40W MED BASE OVEN
39	008231	1	LT FIX CAP AY #VFC-100
40	017631	1	LT FIX SPLICE BOX #VGA-1
41	123001	1	LT FIX VFC-100 MTG PLT/BOX
42	261017	1	LT FIX GLOBE VCGP-100 500F MAX
43	145858	1	LBB1-23-1 DR SEAL FRAME DETAIL
44	145828	1	SHELF PLATED WIRE LBB1-23-1 17.75"X17.75"
45	334111	1	LT FIX 40W G9 OVEN 230V, BJB 77.953.3103.72 W/230V LAMP
46	354826	1	SNAP RING, G9 LT FIXTURE RETROFIT
47	354829	1	FLANGE, HOUSING END, G9 FIXTURE
48	272502	4	SCW MAC PHTX ZIN #08-32X0.50

NOTE: 120V OR 240V MODEL

REV.	DATE	BY	DESCRIPTION	ECO #
D	12/19/2019	JRS	ADDED 7'-3" TO TITLE	
C	6/19/2017	NPL	ADDED ITEMS 47 AND 48. ADDED BOXES AROUND THE 120V OPTION AND THE 240V OPTION. DELETED A 120V NOTE FROM ITEMS 41, 42 AND 44 B/C THEY ARE IN BOTH OPTIONS. CHANGED ITEM 38 FROM PN 008220 TO PN 141894.	1700105
B	9/30/2016	DWW	RE MOVED SWITCH COMPONENTS (ITEMS 24-29). NO SWITCH OPTION FOR 188. ADDED G9 RETROFIT ITEMS 45 (344111) & 46 (354826) & QTY3 OF ITEM 10 (133033).	
A	12/11/14	RWF	INITIAL RELEASE	
			DESCRIPTION	

UNLESS OTHERWISE SPECIFIED DIMENSIONAL TOLERANCES
 DECIMALS .XX ± .030
 .XXX ± .005

THE DESIGN SHOWN ON THIS DRAWING IS THE EXCLUSIVE PROPERTY OF DESPATCH INDUSTRIES, AND IS SUBMITTED WITH THE UNDERSTANDING THAT IT IS NOT TO BE USED FOR ANY PURPOSE OR ITS CONTENTS DISCLOSED WITHOUT THE WRITTEN CONSENT OF DESPATCH INDUSTRIES.

COPYRIGHT ©2010 DESPATCH INDUSTRIES.

DR'N	RWF	SCALE	2:13
CK'D	DLM	DATE	12/11/2014
APP	DLM	SMALL TO	---

PARENT NO.	SERIAL NO.	CUSTOMER NO.
LBB1-23-2/-3 REPLACEMENT PARTS		
D	---	340783
D	---	1 of 1
SIZE	CODE	PART NUMBER
		REV
		SHT OF SHTS