



| LIST OF MATERIAL | | | |
|------------------|-------------|------|---|
| ITEM | PART NUMBER | QTY. | DESCRIPTION |
| 1 | 338063 | 1 | HEATER COVER SHT |
| 2 | 338064 | 1 | DUCT SHT LFT |
| 3 | 338065 | 1 | DUCT SHT RHT |
| 4 | 007778 | 1 | HTR TUTCO 3.6KW 208-240V |
| 5 | 336928 | 1 | SLV 2.50 LG ALUM |
| 6 | 336929 | 1 | MTR MTG LAB LINE .04HP |
| 7 | 133033 | 26 | SCW SM PHTX ZIN #06X0.50 |
| 8 | 012777 | 11 | WSHR FLT 18-8 #06 |
| 9 | 016643 | 1 | CLIP TINNEMAN C17120-017-4 |
| 10 | 145783 | 1 | MTR 00.25HPSODP11523015060 |
| 11 | 015876 | 1 | FAN WH 05.75X01.00X0.500 ALZ CCW |
| 12 | 336959 | 2 | SHIM STRIKE 0.50X0.50 SST |
| 13 | 008199 | 2 | LATCH #350F2C |
| 14 | 272516 | 4 | SCW MAC THTX18-8 #08-32X1.00* |
| 15 | 148048 | 2 | DOOR PULL-AUSTIN #AE1-30357 |
| 16 | 143614 | 4 | SCW TT PHTX ZIN #08-32X0.50 |
| 17 | 145929 | 14 | SCW TT PHTX ZIN 0.25-20X1.00 |
| 18 | 015229 | 1 | ENCL BLWR UF12A23 STH 230V |
| 19 | 018967 | 4 | SCW MAC RHSL ZIN #06-32X2.00 |
| 20 | 008484 | 4 | NUT HX ZIN #06-32 |
| 21 | 012814 | 4 | WSHR LK ZIN #06 SHKPRF INT |
| 22 | 144755 | 1 | T/C JJ38005.5-303.5-25048-006 |
| 23 | | | |
| 24 | | | |
| 25 | | | |
| 26 | | | |
| 27 | | | |
| 28 | | | |
| 29 | 057939 | 1 | STOPPER SLCN #10.5 HI-TEMP* |
| 30 | 132623 | 6 | RIV POP SST .12DIA X.032X.187 LG |
| 31 | 336950 | 1 | PORT TUBE 2.0 IN (6.21LG) |
| 32 | 336951 | 1 | 2 INCH PORT COVER PLATE |
| 33 | 163101 | 1 | T/C K48U-006-00-19-F1B084-6 |
| 34 | 163023 | 1 | POCKET PROBE READOUT |
| 35 | 272527 | 2 | SCW TT PHTX ZIN #08-32X1.00 |
| 36 | 325729 | 1 | METER DG OMGA HH801B HAND HELD DUAL CHANNEL |
| 37 | | | |
| 38 | | | |
| 39 | 017631 | 1 | LIGHT FIXTURE SPLICE BOX #VGA-1 |
| 40 | 123001 | 1 | LT FIX VFC-100 MTG PLY/BOX |
| 41 | | | |
| 42 | 272579 | 4 | SCW MAC THTX18-8 0.25-20X0.50 |
| 43 | 145773 | 1 | LBB2-12-1 DR SEAL FRAME DETAIL |
| 44 | 337057 | 1 | SHELF WIRE PLTD 19.75" X 29.75" LBB2-12 |
| 45 | 334111 | 1 | LT FIX 40W G9 OVEN 230V, BJB 77.953.3103.72 W/230V LAMP |
| 46 | 354826 | 1 | SNAP RING, G9 LT FIXTURE RETROFIT |
| 47 | 354829 | 1 | FLANGE, HOUSING END, G9 FIXTURE |
| 48 | 272502 | 4 | SCW MAC PHTX ZIN #08-32X0.50 |

NOTE: 240V MODEL

| REV. | DATE | BY | DESCRIPTION | ECO # |
|------|------------|-----|---|---------|
| D | 12/19/2019 | JRS | ADDED "7-3" TO TITLE | |
| C | 6/19/2017 | NPL | REMOVED ITEMS 37(PN 008220), 38(PN 008231), AND 41(PN 261017) THAT WERE PART OF THE 120V OPTION. REMOVED LABELING ON ITEMS 45 AND 46 B/C THERE IS ONLY 240V OPTION. ADDED ITEMS 47 AND 48 4X OF 46. | 1700105 |
| B | 10/3/2016 | DWW | REMOVED SWITCH COMPONENTS (ITEMS 23-28). NO SWITCH OPTION FOR LBB. REPLACED ITEM 44 (145779) W/ 337057. ADDED G9 RETROFIT ITEMS 45 (334111) & 46 (354826) & 3X OF ITEM 7 (133033) | UPDATED |
| B | 12/8/14 | RWF | INITIAL RELEASE | |

UNLESS OTHERWISE SPECIFIED
ALL DIMENSIONS IN INCHES:
DIMENSIONAL TOLERANCES:
.XX ± .03
.XXX ± .005
ANGLES ± 1°

THE DESIGN SHOWN ON THIS DRAWING IS THE EXCLUSIVE PROPERTY OF DESPATCH INDUSTRIES, AND IS SUBMITTED WITH THE UNDERSTANDING THAT IT IS NOT TO BE USED FOR ANY PURPOSE OR ITS CONTENTS

COPYRIGHT ©2010 DESPATCH INDUSTRIES

| | | | |
|------|------|------------|-----------|
| DRN | RWF | SCALE | NONE |
| CK'D | DLM | DATE | 12/8/2014 |
| APP | -DLM | SIMILAR TO | SWIS |

Despatch
MINNEAPOLIS, MN • LAKEVILLE, MN

LBB2-12-2/3 REPLACEMENT PARTS

D — 340781 D 1 of 1

SIZE: PART NUMBER REV SHT OF SHTS

©DESPATCH, revised 20100901