



| LIST OF MATERIAL | | | |
|------------------|-------------|------|---|
| ITEM | PART NUMBER | QTY. | DESCRIPTION |
| 1 | 336924 | 1 | HEATER COVER SHT |
| 2 | 336925 | 1 | DUCT SHT LFT |
| 3 | 336926 | 1 | DUCT SHT RHT |
| 4 | 336928 | 1 | SLV 2.50 LG ALUM |
| 5 | 336929 | 1 | MTR MTG LAB LINE .04HP |
| 6 | 133033 | 40 | SCW SM PHTX ZIN #06X0.50 |
| 7 | 012777 | 11 | WSHR FLT 18-8 #06 |
| 8 | 143614 | 4 | SCW TT PHTX ZIN #08-32X0.50 |
| 9 | 016643 | 1 | CLIP TINNERMAN C17120-017-4 |
| 10 | 145929 | 14 | SCW TT PHTX ZIN 0.25-20X1.00 |
| 11 | 145783 | 1 | MTR 00.25HPSODP11523015060 |
| 12 | 015876 | 1 | FAN WH 05.75X01.00X0.500 ALZ CCW |
| 13 | 007818 | 3 | HTR TUTCO 1.2KW 120/240V |
| 14 | 015229 | 1 | ENCL BLWR UF12A23 STH 230V |
| 15 | 018967 | 4 | SCW MAC RHSL ZIN #06-32X2.00 |
| 16 | 008484 | 4 | NUT HX ZIN #06-32 |
| 17 | 012814 | 4 | WSHR LK ZIN #06 SHKPRF INT |
| 18 | 336959 | 2 | SHIM STRIKE 0.50X0.50 SST |
| 19 | 008199 | 2 | LATCH #350F2C |
| 20 | 272579 | 4 | SCW MAC THTX18-8 0.25-20X0.50 |
| 21 | 272516 | 4 | SCW MAC THTX18-8 #08-32X1.00* |
| 22 | 057939 | 1 | STOPPER SLCN #10.5 HI-TEMP* |
| 23 | 132623 | 6 | RIV POP SST .12DIAx.032X.187 LG |
| 24 | 336950 | 1 | PORT TUBE 2.0 IN (6.21LG) |
| 25 | 336951 | 1 | 2 INCH PORT COVER PLATE |
| 26 | 163101 | 1 | T/C K48U-006-00-19-F1B084-6 |
| 27 | 163023 | 1 | POCKET PROBE READOUT |
| 28 | 272527 | 2 | SCW TT PHTX ZIN #08-32X1.00 |
| 29 | 325729 | 1 | METER DG OMGA HH801B HAND HELD DUAL CHANNEL |
| 30 | | | |
| 31 | | | |
| 32 | | | |
| 33 | | | |
| 34 | | | |
| 35 | | | |
| 36 | | | |
| 37 | | | |
| 38 | 017631 | 1 | LT FIX SPLICE BOX #VGA-1 |
| 39 | 123001 | 1 | LT FIX VFC-100 MTG PLT/BOX |
| 40 | | | |
| 41 | 144755 | 1 | T/C JJ38005.5-303.5-25048-006 |
| 42 | 148048 | 2 | DOOR PULL-AUSTIN #AE1-30357 |
| 43 | 145814 | 1 | DR SEAL FRAME DETAIL LBB2-18/27-1 39.75" x 39.50" |
| 44 | 336978 | 2 | 36.75 X 23.75 SHELF PLATED WIRE |
| 45 | | | |
| 46 | 334111 | 1 | LT FIX 40W G9 OVEN 230V, BJB 77.953.3103.72 W/230V LAMP |
| 47 | 354826 | 1 | SNAP RING, G9 LT FIXTURE RETROFIT |
| 48 | 354829 | 1 | FLANGE, HOUSING END, G9 FIXTURE |
| 49 | 272502 | 4 | SCW MAC PHTX ZIN #08-32X0.50 |

NOTE: 240V MODEL

| REV. | DATE | BY | DESCRIPTION | ECO # |
|------|------------|-----|---|---------|
| D | 12/19/2019 | JRS | ADDED "1-3" TO TITLE | |
| C | 6/19/2017 | NPL | REMOVED ITEMS 36 (PN 008220), 37 (PN 008231), AND 40 (PN 261017) THAT WERE PART OF THE 120V OPTION. REMOVED LABELING ON ITEMS 46 AND 47 B/C THERE IS ONLY 240V OPTION. ADDED ITEMS 48 AND 49. | 1700105 |
| B | 9/30/2016 | DWW | REMOVED SWITCH COMPONENTS (ITEMS 30-35) SWITCH NOT LBB OPTION. REPLACED ITEMS 44 (105498) & 45 (016471) W/ 336978 (QTY 2). ADDED G9 RETROFIT ITEMS 46 & 47 & QTY 3 OF ITEM 6. | UPDATED |
| A | 12/12/14 | RWF | INITIAL RELEASE | |

UNLESS OTHERWISE SPECIFIED
ALL DIMENSIONS IN INCHES:
DIMENSIONAL TOLERANCES:
.XX ± .03
.XXX ± .005
ANGLES ± 1°

THE DESIGN SHOWN ON THIS DRAWING IS THE EXCLUSIVE PROPERTY OF DESPATCH INDUSTRIES, AND IS SUBMITTED WITH THE UNDERSTANDING THAT IT IS NOT TO BE USED FOR ANY PURPOSE OR ITS CONTENTS

| | | | |
|------|------|------------|------------|
| DRN | RWF | SCALE | NONE |
| CK'D | DLM | DATE | 12/12/2014 |
| APP | -DLM | SIMILAR TO | SWIS |

Despatch
MINNEAPOLIS, MN • LAKEVILLE, MN

LBB2-18-2/3 REPLACEMENT PARTS

COPYRIGHT ©2010 DESPATCH INDUSTRIES

D 340780 D 1 of 1

SIZE PART NUMBER REV SHIT OF SHITS

c:\despatch\d_revised\20100901