



LIST OF MATERIAL			
ITEM	PART NUMBER	QTY.	DESCRIPTION
1	336970	1	HEATER COVER SHT
2	336971	1	DUCT SHT LFT
3	336972	1	DUCT SHT RHT
4	007819	2	HTR TUTCO 2.4KW 120/240V
5	015876	1	FAN WH 05.75X01.00X0.500 ALZ CCW
6	336928	1	SLV 2.50 LG ALUM
7	016643	1	CLIP TINNEMAN C17120-017-4
8	336929	1	MTR MTG LAB LINE .04HP
9	145783	1	MTR 00.25HPSODP11523015060
10	133033	30	SCW SM PHTX ZIN #06X0.50
11	012777	11	WSHR FLT 18-8 #06
12	145929	14	SCW TT PHTX ZIN 0.25-20X1.00
13	336959	2	SHIM STRIKE 0.50X0.50 SST
14	008199	2	LATCH #350F2C
15	272516	4	SCW MAC THTX18-8 #08-32X1.00*
16	272579	4	SCW MAC THTX18-8 0.25-20X0.50
17	015229	1	ENCL BLWR UF12A23 STH 230V
18	018967	4	SCW MAC RHSL ZIN #06-32X2.00
19	008484	4	NUT HX ZIN #06-32
20	012814	4	WSHR LK ZIN #06 SHKPRF INT
21	144755	1	T/C JJ38005.5-303.5-25048-006
22	057939	1	STOPPER SLCN #10.5 HI-TEMP*
23	132623	6	RIV POP SST .12DIA X.032X.187 LG
24	336950	1	PORT TUBE 2.0 IN (6.21LG)
25	336951	1	2 INCH PORT COVER PLATE
26	163023	1	POCKET PROBE READOUT
27	163101	1	T/C K48U-006-00-19-F1B084-6
28	272527	2	SCW TT PHTX ZIN #08-32X1.00
29	325729	1	METER DG OMGA HH801B HAND HELD DUAL CHANNEL
30			
31			
32			
33			
34			
35			
36			
37			
38	017631	1	LT FIX SPLICE BOX #VGA-1
39	123001	1	LT FIX VFC-100 MTG PLT/BOX
40			
41	148048	2	DOOR PULL-AUSTIN #AE1-30357
42	143614	4	SCW TT PHTX ZIN #08-32X0.50
43	145814	1	DR SEAL FRAME DETAIL LBB2-18/27-1 39.75" x 39.50"
44	337058	1	36.75 X 36.75 SHELF PLATED WIRE
45	334111	1	LT FIX 40W G9 OVEN 230V, BJB 77.953.3103.72 W/230V LAMP
46	354826	1	SNAP RING, G9 LT FIXTURE RETROFIT
47	354829	1	FLANGE, HOUSING END, G9 FIXTURE
48	272502	4	SCW MAC PHTX ZIN #08-32X0.50

NOTE: 240V MODEL

REV.	DATE	BY	DESCRIPTION	ECO #
D	12/19/2019	JRS	ADDED "7-3" TO TITLE	
C	6/16/2017	NPL	REMOVED 120V COMPONENTS (ITEMS 36, 37, 40) FROM DRAWING AND BOM. ADDED ITEMS 47 AND 48.	1700105
B	10/3/2016	DWW	REMOVED SWITCH COMPONENTS (ITEMS 30-35). SWITCH NOT LBB OPTION. REPLACED ITEM 44 (016472) W/ 337058. ADDED G9 RETROFIT ITEMS 45 & 46. ADDED QTY 3 OF ITEM 10 (133033).	UPDATED
B	12/12/14	RWF	INITIAL RELEASE	

UNLESS OTHERWISE SPECIFIED DIMENSIONAL TOLERANCES
 DECIMALS .XX ± .030
 .XXX ± .005

THE DESIGN SHOWN ON THIS DRAWING IS THE EXCLUSIVE PROPERTY OF DESPATCH INDUSTRIES. IT IS SUBMITTED WITH THE UNDERSTANDING THAT IT IS NOT TO BE USED FOR ANY PURPOSE OR ITS CONTENTS DISCLOSED WITHOUT THE WRITTEN CONSENT OF DESPATCH INDUSTRIES.

COPYRIGHT ©2014 DESPATCH INDUSTRIES

DR'N	RWF	SCALE	1:10
CK'D	DLM	DATE	12/12/2014
APP	DLM	SMALL TO	---

Despatch
 MINNEAPOLIS, MN • KEVILLEN, MN

LBB2-27-2/-3 REPLACEMENT PARTS

D 340779 D 1 of 1

SIZE	CODE	PART NUMBER	REV	SHT OF SHTS
------	------	-------------	-----	-------------