



		PART MATERIAL		
ITEM	PART NUMBER	QTY.		DESCRIPTION
1	338640	1		LCC1-51-5 SHELF SUPPORT/DUCT WLDMT
2	156429	1		DOOR PULL
3	016008	4		BRG SLVFLG BRZ 0.38X0.25X0.50*
4	185436	1		GAGE MAGHLC 2310-SF 5-0-5
5	098365	1		FAN AY EXH LAB
7	098168	1		HTR TUTCO 6.0KW SST 208/240V
8	156403	1		MTR 00.25 HP TENV208230160 CFACE
9	156426	1		LATCH BRACKET SST LCC
10	156433	1		HNG PIVOT LOWER LEFFT ASSY PLATED
11	156434	2		HNG PIVOT UPPER LEFT ASSY PLATED
12	156435	1		HNG PIVOT UPPER RIGHT ASSY PLATED
13	345266	1		FAN WHEEL 6.25 X 1.50 X .500 SST CCW
14	162716	2		SEAL 2.0 ODX0.44 IDX0.12 GORE-TEX (SOFT)
15	122059	1		LATCH STHCO#A7-10-301-20 DRAW
16	326699	3		VLV CHK 62-502-01 0.38"NP SST
17	085674	1		VLV NDL #N10B .125 BRASS DELTROL
18	051706	1		METER-FLW RMA-10D-BV 20-200SCFH/5-95 LPM WITH VALVE
19	141810	1		VLV SOL 8210G001HW .38 NC 240
20	141812	1		VLV SOL SC8262H202 .25 NC 240
21	156463	1		METER-FLW VFB-86-SSV .5-5 GPM
22	156464	1		VLV SOL 8316G54 .38 NC 240V 3WAY
23	160240	1		VALVE RELIEF VENT LCC1-51N
24	308437	1		COIL WATER COOL 5X12.88-2R-38/168
25	158802	1		LCC1-51 DOOR SEAL FRAME DETAIL-SILICONE
26	156495	1		LCC1-51 MOLDED VITON DOOR SEAL FRAME
27	156488	1		FLTR HEPA #00FH-18Z18Z06-5U-2-E-A-2D-O/OC
28	156457	2		FILTER MTG ANGLE
29	160226	1		FILTER HEPA #4P4 TERMIKFI
30	337739	1		THREAD RD BRASS .38-16X2.5"
31	156482	2		SHELF LCC1-51 SST WIRE 19.75 X 22.62
32	163405	4		NUT HX BRASS #06-32
33	010891	1		SCW MAC RHSL 18-8 #06-32 X 0.75

DR'N	DWW	SCALE	1:8
CK'D	TJK	DATE	6/13/2016
APP	GDY	SMALL TO	---

Despatch
MINNEAPOLIS, MN • LAKEVILLE, MN

LCC1-51-5 REPLACEMENT PARTS

UNLESS OTHERWISE SPECIFIED DIMENSIONAL TOLERANCES

DECIMALS
.XX ± .030
.XXX ± .005

THE DESIGN SHOWN ON THIS DRAWING IS THE EXCLUSIVE PROPERTY OF DESPATCH INDUSTRIES. IT IS SUBMITTED WITH THE UNDERSTANDING THAT IT IS NOT TO BE USED FOR ANY PURPOSE OR ITS CONTENTS DISCLOSED WITHOUT THE WRITTEN CONSENT OF DESPATCH INDUSTRIES.

COPYRIGHT ©2014 DESPATCH INDUSTRIES

A	6/13/2016	DWW	INITIAL RELEASE	
REV.	DATE	BY	DESCRIPTION	ECO #

D	---	---	351992	A	1 of 1
SIZE	CODE	PART NUMBER	REV	SHT OF SHTS	