



- NOTE:
- 1- ALUMINIZED BODY ONLY
 - 2- STAINLESS STEEL BODY ONLY
 - ▷ RAD PARTS ONLY
 - ▷ RFD PARTS ONLY

ITEM	PART NUMBER	QTY.	DESCRIPTION
1	030327	8.0 FT	SEAL TAPE FBRGLS .25X2.00
2	355092	2	DR. CATCH, R_D1-42
3	308714	2	LATCH SLIDE BOLT-MACHINED
4	315413	2	LATCH HSG #5658-50
5	147588	1	LATCH 2PT CTR CASE ZN#5646-50
6	147542	1	HANDLE ASSY DII #3 BLK PWDR COAT
7	160872	2	HINGE ASSY #4DII HIGH HEEL
8	044908	12	SCW TT HWPO STL 0.31-18X1.00
9	272556	6	SCW TT HW ZIN 0.25-20X0.50
10	272526	18	SCW DT PHTX ZIN #08-18X0.63
11	044922	2	SCW TT SKTXPO STL #10-32X.50
12	020702	8.0 FT	SEAL #GI-T15CB55 0.75X2.50KEV *
13	305576	2	DISC SPRING "ROD SIZE 5/16"
14	308713	2	PIN COTTER .094DIA X .75L ZINC
15	044703	18	SCW DT THTX STL .25-20 X 1.50
16	122185	22	SCW TT PHTX ZIN #08-32 X 0.38
17			
18	018386	1	PIPE CAP 304 1.00
19	012797	4	WSHR FLT ZIN 0.25
20	024690	1	SFT SEAL PLT 6.5 SQ 18GA .75D
21	147571	1	FAN RING INL 05.5 FC 304SST
22	147580	1	FAN WH 5.50 X 3.00 X .625 SST CCW
23	010836	4	SCW CAP HH ZIN 0.31-18X1.00
24	012798	4	WSHR FLT ZIN 0.31
25	012837	6	WSHR LK ZIN 0.31 SPG
26	081308	1	SFT SEAL 6.0X6.0X0.12 FBGLS PATCH
27	272514	4	SCW MAC THTX18-8 #08-32X0.50
28	272524	4	SCW DT HW ZIN #08-18X0.75
29	147570	1	MTR 00.50HPSODP115230 561*
30	146549	1	SLV 2.0DIA X .065WALL 4.44 SST R_D1-42
31	013569	1	T/C FIG 1/8NPT FOR 3/16OD T/C
32	015229	1	ENCL BLWR UF12A23 3TH 230V
33	017800	1	T/C JJ38E-018-00-16-F1 A060-024-2
34			
34B	149807	2	PLT EXPL RLF CVR 304SSTRAD1-42
35	008484	4	NUT HX ZIN #06-32
36	012814	4	WSHR LK ZIN #06 SHKPRF INT
37	018967	4	SCW MAC RHSL ZIN #06-32X2.00
38	142861	1	CORD BLWR PWR GDTEC HP50-60DESP [COILED AND TIE WRAPPED]
39	141806	2	TERM BLK PHNX MSTB2.5/3-ST-5.08 24-12 GA 10A PLUG-IN 3 POSITION PHNX 1757022
41	146252	1	HTR AY RAD/RFD1-42
42	012833	8	WSHR LK ZIN #08 SPG
43			
43B	146529	1	PLATE INSIDE FLOOR R_D1-42 (SST)
44			
44B	146550	1	SUPPLY BAFFLE R_D1-42(SST)
45			
45B	146552	1	RETURN BAFFLE R_D1-42(SST)
46			
46B	147575	2	FRAME EXP RLF 4.0 X 19.0 .53SQ RFD1-42 (SST)
47	149803	1	PAN FAN RFD1-42 EXH
48	111535	2	SHT SST 321 .002TH X 24"W X 600 IN*
49	407373	1	DMPR SLIDE 2.50 DIA FA .56 OPN W/ 0.56 HOLE
50	187702	4	SCW CAP PHSL BRS #08-32X0.25
51	032565	4	SCW DT HHSK STL .25-14 X 1.50
52	143614	4	SCW TT PHTX ZIN #08-32X0.50
53	012778	4	WSHR FLT 18-8 #08
54	027452	2	TU FTG .25TX.12PME, EL, PUSH, SWIVEL
55	051281	4	TU FTG B400-2-2 EL .25TX.12PM
56	012544	1	TU FTG BRS COMP .25TX.12P ELL*
57	012545	1	TU FTG#68M-04X02 .25TX.12P ST*
58	258994	2	TU FTG 3014.56.11 .25TX.12PPT *
59	143482	2	SW AF JD-2 .10-4.00 IN WC DIN
60	345667	1	TUBE AFS PRESSURE
61	345668	1	TUBE AFS VACUUM
62	158721	1	TU AFS VACUUM, RFD1-42
63	158720	1	TU AFS PRESSURE, RFD1-42
64	145593	2	CONN DIN MINCTL 9067-128
65	027380	1	TU NYL 1094P56-00 LEGRIS 0.25
66	015787	1	SHELF PLTD WIRE 17.75 X 19.62 (R_D1-42)
67	149804	1	DMPR AY RAD1-42 STD EXH
68	058667	2	PIPE PLUG 304 0.12 HH
69	149809	1	DMPR SLIDE 2.50 DIA FA 0 OPN RAD1-42
70	112667	1	FAN AY FASCO MODIFIED
71	146541	2	ROD DOOR LATCH R_D1-42
72	030322	.17 FT	SEAL TAPE FBRGLS .12 X1.50*
73	030317	0.4 FT	SEAL TAPE FBRGLS .06 X2.00*
74	010834	2	SCW CAP HH ZIN 0.31-18X0.50
75	008485	8	NUT HX ZIN #08-32
76	149798	1	FE RFD1-42 0.025 GPH MEK
77	272563	4	SCW CAP HW ZIN 0.25-20X0.75*
78			
82	272501	6	SCW MAC PHTX ZIN #08-32 X 0.25

REV.	DATE	BY	DESCRIPTION	ECO #
D	9/13/2023	JRS	ITEM 7, 146872 WAS 338703	UPDATE
C	8/11/2022	JRS	REMOVED ITEMS: 17, 272501 QTY 4, ITEM 34A, 402434 QTY 2, ITEM 45A, 402498 QTY 1, ITEM 44A, 402410 QTY 1, ITEM 45A, 402409 QTY 1, ITEM 46A, 402435 QTY 2	UPDATE
B	10/28/2020	BMW	REPLACED FA DMPR SLIDE	2000141
A	8/2/2019	HEV	INITIAL RELEASE	

UNLESS OTHERWISE SPECIFIED DIMENSIONAL TOLERANCES

DECIMALS .XX ± .030 .XXX ± .005

THE DESIGN SHOWN ON THIS DRAWING IS THE EXCLUSIVE PROPERTY OF DESPATCH INDUSTRIES. IT IS SUBMITTED WITH THE UNDERSTANDING THAT IT IS NOT TO BE USED FOR ANY PURPOSE OR ITS CONTENTS DISCLOSED WITHOUT THE WRITTEN CONSENT OF DESPATCH INDUSTRIES.

COPYRIGHT ©2014 DESPATCH INDUSTRIES

DR'N	HEV	SCALE	1:8
CK'D	MJS	DATE	08/02/2019
APP	MJS	SMALLER TO:	346353

Despatch
MINNEAPOLIS, MN • WAKEFIELD, MA

R_D1-42-XE1 REPLACEMENT PARTS

D	404425	D	1 of 1
SIZE	CODE	PART NUMBER	REV
			SHT OF SHTS